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Ideal solution for dehydrated and damaged hands and feet skin. The skin will be moisturized and smoothed. Paraffin will strengthen the nails, stimulate the blood circulation and work therapeutically for bones and joints.

PACKING LIST:

- Paraffin heater with protective grid on the bottom of device;
- 2x paraffin with a pencil for paraffin;
- 1 pair of gloves (frotte), 1 pair of socks (frotte) and PE gloves 25 pcs.

OPERATION INSTRUCTION:

- Plug in the device and put paraffin in a heater.
- Put switch on HIGH position.
- Paraffin warm up in about 20/30 minutes.
- Then put switch on LOW position. In this mode thermostat will keep appropriate temperature (ideal to make a procedure).
- Before you start the procedure, put a protective grid at the bottom of device.
- When you end the procedure, please put switch on OFF position.

TECHNICAL INSTRUCTION:

- Power: 200W
- Voltage: AC 220-240V ~ 50/60Hz
- Dimensions: 33cm x 24,5cm x 18cm
- Dimensions of vessel: 26cm x 18cm x 9,5cm
- Thermostat

CAUTIONS:

- Please read the instruction carefully before using the device.
- Please use the device according to its destination.
- Please unplug the power supply, when you finish using the device.
- Disconnect the device from the mains during thunderstorms or when it is not to be used for an extended period of time.
- Do not use the device continuously.
- The device could be used by mentally, physically disabled people (also children) or by people without appropriate experience or knowledge only if suitable supervision or instruction will be ensured. Keep away from children. Children without the adult's supervision should not clean and maintain the device.
- The inseparable power cord should be changed only in manufacturer's company or by an employee of the service or by a qualified person.

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